



# IndiMag 48s

Fast and versatile nucleic acid extraction

- Automated, magnetic bead-based extraction for 1 to 48 samples
- Just 30 minutes per run
- Open software for easy protocol customization
- Effective decontamination using UV light
- Compact footprint to fit any lab
- Optimized plastic use thanks to 3 flexible cartridge formats



## IndiMag Cartridges

Prefilled for maximum efficiency

- Save significant hands-on time and plastics with IndiMag prefilled reagent cartridges
- Simply add the lysis buffer and your samples
- Available in 8 and 24 sample sizes for IndiMag 48/48s
- Proven high-throughput performance

**Get in touch for more information**

Email us: [support@indical.com](mailto:support@indical.com) | Call us: **+49 341 124 54 0** | Learn more at: [www.indical.com](http://www.indical.com)

## IndiMag 48s and compatible plasticware

Product	Short description	Cat. no. (number of preps)
IndiMag 48s (100~240V)	Automated extraction system for rapid isolation of nucleic acids using magnetic bead technology	IN943048s
IndiMag 48s Service Solutions	Extended warranty for 2–4 years for maximum productivity and instrument uptime <small>(Please contact us for further details)</small>	IMG48SERV
IndiMag 48 PW 1-Sample Block	Plasticware (4-well blocks) for 368 x 1 samples on IndiMag 48/48s	PW940123 (368)
IndiMag 48 PW 8-Sample Blocks	Plasticware (32-well blocks) for 70 x 8 samples on IndiMag 48/48s	PW940166 (560)
IndiMag 48 PW 24-Sample Blocks	Plasticware (96-well blocks) for 28 x 24 samples on IndiMag 48/48s	PW940187 (672)
IndiMag 48 PW Rod Cover	Rod cover (8-strips) for 84 x 8 samples on IndiMag 48/48s	PW940237 (672)

## IndiMag 48 and IndiMag 48s spare parts

Product	Cat. no. (pcs)		Product	Cat. no. (pcs)
	For IndiMag 48	For IndiMag 48s		IndiMag 48s
UV lamp	1300001 (1)	1300002 (1)	Cartridge adapter	1300015 (2)
Fuse	1300003 (2)	1300004 (2)	Power cord US	1300016 (1)
Magnet rod	1300005 (1)	1300006 (1)	Power cord CN	1300017 (1)
Calibration plates	1300007 (2)	1300008 (2)	Power cord UK	1300018 (1)
Filling adapter	1300009 (1)	1300010 (1)		
Adapter	1300011 (1)	1300012 (1)		
Magnet stand	1300013 (1)	1300014 (1)		

## IndiMag nucleic extraction kits and prefilled reagent cartridges

Product	Short description	Cat. no. (number of preps)		
		S	M	L
NEW! IndiMag Pathogen IM48 Cartridge	Prefilled reagent cartridges for automated extraction of nucleic acids with IndiMag 48 and IndiMag 48s	SP947654P608 (6 x 8)	SP947654P224 (2 x 24)	
IndiMag Pathogen Kit w/o plastics	Automated extraction of nucleic acids with IndiMag 48 and IndiMag 48s <small>(Please order plasticware separately)</small>			SP947257 (384)
IndiMag Mastitis Kit	Automated extraction of bacterial DNA from milk samples with IndiMag 48 and IndiMag 48s <small>(Please order plasticware separately)</small>			SP947757 (384)

Learn more at: [www.indical.com](http://www.indical.com)

Email us: [support@indical.com](mailto:support@indical.com) | Call us: +49 341 124 54 0 | Fax us: +49 341 124 54 60 | Follow us: [on LinkedIn](#)

For up-to-date licensing information and product-specific disclaimers, see the respective handbook or user manual. Regulatory requirements vary by country, products may not be available in your geographic area. Trademarks: INDICAL®, IndiMag® (INDICAL BIOSCIENCE GmbH). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-IMG48s-IM48-EN-20210428 © 2021 INDICAL BIOSCIENCE

**INDICAL**  
BIOSCIENCE

INDICAL BIOSCIENCE GmbH  
Deutscher Platz 5b  
04103 Leipzig, Germany

